



Z-Band UHF RFID Wristbands

Improve Accuracy and Operational Efficiency

Zebra's Z-Band® UHF RFID durable direct thermal printable wristband solutions simplify the identification of patients in clinical environments or guests at event venues. These disposable RFID wristbands are specifically designed for optimal use with Zebra's ZD500R desktop printer and Zebra RFID readers. Zebra's UHF RFID wristbands offer a wristbanding solution to fit your RFID application, helping you better manage events or ensure positive patient identification.

Z-Band Direct RFID SR Wristband

Zebra's Z-Band Direct RFID SR wristband is best for use cases where the wristband will be in close proximity to the RFID reader.

Common Applications

Healthcare*

- At the patient bedside through a blanket
- During a patient's surgery through the drape
- Allows caregivers to obtain immediate access to critical patient records and update them at the bedside
- Quickly associate supplies and equipment used in treatments

Hospitality

- Guest Identification
- Access Management
- Cashless POS



What Can Passive UHF RFID Do?

Passive UHF RFID technology can utilize short and long range RFID readers to offer the tracking of any tagged asset with last known read location throughout a facility.

Part Number	Color	Size	Inlay	Read Range	Latex-free
3014578	White	1" x 11"	BT0600	Up to 3m	Yes

FACT SHEET

Z-BAND UHF RFID WRISTBAND

Z-Band UltraSoft RFID LR Wristband

Zebra's Z-Band UltraSoft LR Wristband is uniquely designed with a specialty inlay designed for use against the human body. The design allows the wristband to flag away from the body, enabling longer read ranges. Ideal for use when overhead readers or portals are being used to provide best last known location tracking abilities.



Common Applications

Healthcare*

- Streamlining patient tracking throughout a medical facility

Hospitality

- Guest tracking throughout a venue

Part Number	Color	Size	Inlay	Read Range	Latex-free
ZIPRD3015155	White	1.1875" x 14.25"	ZBR2002	Up to 17m **	Yes

Custom Solutions

Zebra can create a custom wristband featuring your preferred inlay and chip combination, as well as pre-printed logos and colors.

*Wristband should be removed from the patient during an MRI.

**Read range is highly contingent on where inlay is placed in relation to the body.

For more information, visit www.zebra.com/wristbands



NA and Corporate Headquarters
+1 800 423 0442
inquiry4@zebra.com

Asia-Pacific Headquarters
+65 6858 0722
contact.apac@zebra.com

EMEA Headquarters
zebra.com/locations
contact.emea@zebra.com

Latin America Headquarters
+1 786 245 3934
la.contactme@zebra.com